

| Sr. No. | Item Name | Pak Size | Qty. |
|---------|---|--------------|------|
| 1 | Micro spin magnetic stirring bar 2 into 6mm, 10 | Pack of 10 | 2 |
| 2 | Octagon Magnetic Stirrer Bar,PTFE 6 into 10mm, 10 | Pack of 10 | 5 |
| 3 | Octagon Magnetic Stirrer Bar,PTFE 8 into 14mm, 10 | Pack of 10 | 5 |
| 4 | Octagon Magnetic Stirrer Bar PTFE 8 into 22mm, 10 | Pack of 10 | 5 |
| 5 | Racked graduated tips 200µL, 1000 | Pack of 1000 | 10 |
| 6 | Micro Tips Bulk 10µL, 1000 | Pack of 1000 | 10 |
| 7 | Micro Centrifuge Tube 5mL, 100 | Pack of 100 | 10 |
| 8 | Micro Centrifuge Tube 2mL, 500 | Pack of 500 | 5 |
| 9 | Micro Centrifuge Tube zero point 5mL, 1000 | Pack of 1000 | 2 |
| 10 | Dessicators 300 mm, 1 | 01 pack | 2 |
| 11 | Pipette rack horizontal, 6 | Pack of 06 | 2 |
| 12 | Pipettor stand, 1 | 01 Pack | 1 |
| 13 | Sample container 57 into 73 100ml, 280 | Pack of 280 | 1 |
| 14 | Self standing centrifuge tube 50ml, 500 | Pack of 500 | 2 |
| 15 | Tip refill system for 10 µl, 10 | Pack of 10 | 1 |
| 19 | PIPET-Variable volume pipette 1000- 10000µL, 1 | 01 pack | 2 |
| 17 | Macro tips 10 ML, 100 | Pack of 100 | 5 |
| 18 | Pasteur pipette 1 ml, 500 | Pack 500 | 10 |

| Sr. No. | Size | Item Description | Specification |
|---------|-----------|---|--|
| 1 | 2*6mm | Microspin magnetic stirring bar 2*6mm | PTFE/ALNICO V Magnet |
| 2 | 6*10mm | Octagon Magnetic Stirrer Bar,PTFE 6*10mm | PTFE/ALNICO V Magnet |
| 3 | 8*14mm | Octagon Magnetic Stirrer Bar,PTFE 8*14mm | PTFE/ALNICO V Magnet |
| 4 | 8*22mm | Octagon Magnetic Stirrer Bar,PTFE 8*22mm | PTFE/ALNICO V Magnet |
| 5 | 1000/Pack | Racked graduated tips 200µL | Virgin hydrophobic polypropylene, conforming to US FDA21 CFR, Autoclavable, Rnase, Dnase and pyrogen free, free of lubricant, dyes and heavy metals. |
| 6 | 1000/Pack | Micro Tips Bulk 10µL | Virgin hydrophobic polypropylene, conforming to US FDA21 CFR, Autoclavable, Rnase, Dnase and pyrogen free, free of lubricant, dyes and heavy metals. |
| 7 | 100/Pack | Micro CentrifugeTube5mL | PP,confirming to USP class VI, Free of datable pyrogens, Dnase, and Rnase. Provide good clarity for viewing content with stand centrifugation upto 20000xg |
| 8 | 500/Pack | Micro CentrifugeTube2mL | PP,confirming to USP class VI, Free of datable pyrogens, Dnase, and Rnase. Provide good clarity for viewing content with stand centrifugation upto 20000xg |
| 9 | 1000/Pack | Micro CentrifugeTube0.5mL | PP,confirming to USP class VI, Free of datable pyrogens, Dnase, and Rnase. Provide good clarity for viewing content with stand centrifugation upto 20000xg |

| | | | |
|----|----------|---|---|
| 10 | 300mm | Dessicators | Thick walled vaccum desicator made of premium plastic, light in weight,PP/PC |
| 11 | 1 | Pipette rack horizontal | 12 PLACES, horizontal, PP, Autoclavable |
| 12 | Place 5 | Pipettor stand | Made of PMMA material light weight, transparent |
| 13 | 280/pack | Sample container 57x73 [100ml] | PP with HDPE closure, ultra clear and unbreakable, leak proof screw cap material confirm to US FDA21 CFR. |
| 14 | 500/pack | Self standing centrifuge tube [50ml] | PP/HDPE Confirming to usp class vi Autoclavable, free of detectable pyrogen, Rnase, Dnase, MaxRCF 3500Xg |
| 15 | 10 µL | Tip refill system for 10 µl | Graduated,PP CONFIRMG TO US FDA21CFR, Autoclavable, free of detectable Dnase Rnase and pyrogen. |
| 16 | 1 | PIPET-Variable volume pipette 1000-10000µL | Rugged construction, made of PVDE/ABS/SS,Smooth piston, individual calibration certificate |
| 17 | 100/Pack | Macrotips 10 ML | PP.Free of detectable D nase m Rnas and pyrogen. |
| 18 | 1ml | Pasteur pipette 1 ml | Chemical Resistance, Graduated, Made of LDPE, Good Clarity for easy liquid level, Suitable for drawing dispensing, transforming and mixing various volumes of liquid. |

* Supplier should produce authorization certificate from their principal company.

Note: Bidders shall read the Specifications/Terms and Conditions carefully and supply the complete consignment on or before the due date. In case of part supply and/or supply after due date LD charges shall be applicable as per GeM Terms & Conditions